

1986

ANNUAL QUALITY ASSURANCE PERFORMANCE REPORT

SECTION 5

LANDFILL LEACH SAMPLES

INORGANIC TRACE CONTAMINANTS SECTION

MOE  
ANN  
ITC  
ALRN

c.1  
a aa

*ALRN*  
G. C. RONAN, DIRECTOR  
Laboratory Services Branch  
Ministry of the Environment

### Copyright Provisions and Restrictions on Copying:

This Ontario Ministry of the Environment work is protected by Crown copyright (unless otherwise indicated), which is held by the Queen's Printer for Ontario. It may be reproduced for non-commercial purposes if credit is given and Crown copyright is acknowledged.

It may not be reproduced, in all or in part, for any commercial purpose except under a licence from the Queen's Printer for Ontario.

For information on reproducing Government of Ontario works, please contact ServiceOntario Publications at [copyright@ontario.ca](mailto:copyright@ontario.ca)

1986

ANNUAL QUALITY ASSURANCE PERFORMANCE REPORT

SECTION 5

LANDFILL LEACH SAMPLES

INORGANIC TRACE CONTAMINANTS SECTION

D G STURGIS  
and J C HIPFNER (editors)

Inorganic Trace Contaminants Section  
Laboratory Services Branch  
Ministry of the Environment

© 1989 HER MAJESTY THE QUEEN IN RIGHT OF ONTARIO

ALRN

INORGANIC TRACE CONTAMINANTS SECTION

SUMMARY

I. Introduction

The Inorganic Trace Contaminants Section of the Ministry of the Environment, Laboratory Services Branch is responsible for the analysis of a wide variety of sample types for metals and non-metals. The use of sensitive instrumentation and methodologies appropriate to the sample matrix, combined with quality assurance programs, ensures that the Section is able to maintain a high standard of analytical performance. This performance is monitored through regular internal quality control and assurance programs as well as participation in interlaboratory round-robins. This QA report summarizes the methodologies used for analysis of these samples and the supporting internal quality assurance data.

This report is assembled in sections that reflect the analyses performed on different sample matrices in support of the programs of the Ministry of the Environment. Coincidentally, these divisions also reflect the supervisory responsibilities within the Section.

II. Quality Control and Assurance

The objectives of the quality control and assurance programs are to ensure that all of the components of the analytical process are under control and to ensure immediate detection and correction of unacceptable analytical performance. The program monitors all of the reagents, instrumentation, calibration and recovery components of the analytical system.

A. Quality Control

Quality control of the analytical process takes place at the instrument level and is intended to ensure that the instrumentation is operating according to established criteria. This control function ensures that instrument calibration, standardization, slope and intercept, and instrumental drift meet these criteria.

B. Quality Assurance

Quality assurance of the analytical process takes place after the results have been generated and is intended to ensure that the analytical protocols of sample preparation and digestion have been carried out correctly. This control function ensures that reagent blanks, digested standards, sample duplicates and recovery materials meet established response criteria.

### III. Report Format

The report consists of one page method summaries and one page data summaries of blanks, between-run controls and within-run duplicates in formats that are common to all of the parameter/matrix combinations. The method summaries give a brief outline of the sample preparation and measurement procedures. The data summaries consist of annual mean values with standard deviations.

For the within-run duplicates, the data set is subdivided into ranges approximating 0 to 20 %, 20 to 50 % and 50 to 100% of the analytical range. All results for duplicates reported to the data base that are "<" or that have been diluted into the range are excluded from the statistical analysis.

The standard deviations for blanks and between-run controls are calculated using formula I. Formula II is used for the calculations for within-run duplicates.

$$sd = \sqrt{[(\sum x^2 - (\sum x)^2)/n/(n-1)]} \dots\dots I$$

$$sd = \sqrt{(\sum d^2/2n)} \dots\dots II$$

where : x = the individual values; n = the number of events  
d = the differences between pairs of duplicates

The data is stored in a personal computer using BMB Manager II files. All data manipulations, reports generated etc, are performed using applications written in Manager Math.

## TABLE OF CONTENTS

		PAGE NUMBER
SUMMARY	INTRODUCTION	i
	QUALITY CONTROL AND ASSURANCE	i
	REPORT FORMAT	ii
	TABLE OF CONTENTS	iii
SECTION 1	APIOs SAMPLES	1.1
	Lovol Filters	1.3
	Precipitation	1.25
	Precipitation Bag leach	1.59
SECTION 2	BIOMATERIAL SAMPLES	2.1
SECTION 3	AIR EMISSION SAMPLES	3.1
	Hivol Filters	3.5
	Dichotomous Filters	3.45
	Dustfall Jars	3.95
SECTION 4	LIQUID INDUSTRIAL WASTE SAMPLES	4.1
SECTION 5	LANDFILL LEACH SAMPLES	5.1
SECTION 6	SOIL LEACHATE SAMPLES	6.1
SECTION 7	DRINKING AND SURFACE WATER SAMPLES	7.1
SECTION 8	SEDIMENT AND SOIL SAMPLES	8.1
	Sediment	8.3
	Soil	8.73
SECTION 9	MUNICIPAL WASTE SAMPLES	9.1
	Raw Sewage	9.3
	Final Effluent	9.35
	Sludge	9.67
SECTION 10	VEGETATION SAMPLES	10.1

## 5. Landfill

### 5.1 Landfill Leach Samples

Landfill samples are samples collected in bore holes and surface runoff in and around landfill sites in order to monitor the materials that may be leaching from the fill site. In contrast to soil leachate samples, these samples are a result of natural leaching of landfill and other material.

Table 5.1 summarizes the parameters determined, the preparation methods used and the instrument types used for the analysis of landfill samples.

TABLE 5.1

Parameter	Collection Device	Preparation	Analysis
Metals	Plastic or glass jars	Acid digest	ICP - AES
Anions	Plastic or glass jars	Filtration	IC
Cyanide	Plastic or glass jars	Distillation	Colorimetry
Sulfide	Plastic or glass jars		Colorimetry
Mercury	Plastic or glass jars	Acid digest	Cold vapour AAS
Hydride	Plastic or glass jars	Acid digest	Hydride AAS

### 5.2 Landfill Quality Assurance

Sub aliquots of landfill samples are analysed separately to generate duplicate results. Blanks consist of the digestion acid or distilled water as appropriate.

Landfill samples frequently have solid, aqueous and organic liquid phases. In general only the aqueous phase is sampled for inorganic analyses, and the QA materials reflect this position. Composite QA samples are normally spiked to achieve measureable levels of all components determined.

Table 5.2 summarizes the QA materials used when landfill samples are being analysed.

TABLE 5.2

Sample Designation	Type	Parameter
qcal,qca2	Standard solution	Cyanide
qcd	Standard solution (0.2 mg/L)	Cyanide
QCLF-3	composite landfill (prep 1985)	ICP Metals
475-3	EPA solution	Hydride Metals



ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Aluminum      TEST CODE: ALUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-100%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 50 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.08 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.03 - 5.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A

mean	2.30 mg/L
std. dev.	0.62 mg/L
R.S.D.	27.2 %

B

Precision of Duplicates-	low range	mid range	high range
s.d.	0.035	0.793	0.036
mean	0.251	1.549	3.577

W 0.05 mg/L

T 0.25 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - assumed 100% as material originally in solution.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

ALUMINUM

IN LANDFILL

Operating Range = 0.030 to 5.0 mg/L

## IN - RUN DUPLICATES

Range	<0.030	0.030 to1.00	1.00 to2.50	2.50 to5.0	>5.0
no.	34	86	10	3	4
s.w.		0.0348	0.7929	0.0357	
mean		0.2511	1.5494	3.5769	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	2.301	0.6249	27.15

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	116	0.108	0.0711

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: ARSENIC  
UNIT: Biomaterials

TEST CODE: ASUT

SAMPLE TYPE: Landfill Lch.

SUPERVISOR: R. Sadana

METHOD CODE: 510BF3

REVISION NO:

DATE: January, 1983

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 100 ml

Container- Glass bottle with bakelite screw cap (16 oz.)

Preservative- 1ml conc. HNO<sub>3</sub> for sample filling 16 oz. bottle

Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-yes % Extracted-

Procedure- Digest samples by mixing 20 ml of sample and 4 ml of acid 6:3:1 (HNO<sub>3</sub>:HClO<sub>4</sub>:H<sub>2</sub>SO<sub>4</sub>) and heating on a hot plate until dense white fumes are formed. Allow sample to cool and add .5 ml of H<sub>2</sub>O and 2.5 ml HCl. Transfer the contents of the beaker to a test tube graduated at the 20 ml mark and bring to volume with washing. Shake and analyze by automated hydride-F.A.A.S. technique.

INTERFERENCES: Excessive concentrations of Cu, Fe, Ni

REPORTING RESULTS: mg/L-2 dec. if <10, 1 dec. if 10-100, 0 dec. if >100

INSTRUMENTATION: Atomic Absorption Spectrophotometer (Varian 1200) with strip chart recorder, peristaltic pump, auto-sampler, open-ended and heated quartz 'T' cell (0.6x10 cm), and gas-liquid separator.

Calibration Range: 0 to 40 ng/ml

Resolution: 0.01 absorbance

Sensitivity: 20 ng/ml gives 0.150 Abs.

Instrument Detection Limit: 1 ng/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.001 to 0.04 mg/L

Accuracy-

Precision of Controls-

	A	B
mean	.412 mg/L	
std. dev.	.015 mg/L	
R.S.D.	3.7 %	

Precision of Duplicates-	low range	mid range	high range
s.d.	0.0002	0.0007	
mean	0.002	0.009	0.669

W .001 mg/L

T .005 mg/L

CONTROL LIMITS:

REMARKS:

# SUMMARY REPORT OF QUALITY CONTROL DATA

ARSENIC

IN LANDFILL

Operating Range = 0.001 to 0.04 mg/L

## IN - RUN DUPLICATES

Range	<0.001	0.001 to 0.01	0.01 to 0.02	0.02 to 0.04	>0.04
no.	18	19	1	1	2
s.w.		0.0002	0.0007	0.0000	
mean		0.0020	0.0090	0.0320	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
475-3	41	0.412	0.0154	3.74

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
BLK	0	0.000	0.0000

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Barium                      TEST CODE: BAUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AAO  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-100%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.008 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.005 - 2.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

	A	B
mean	.215 mg/L	
std. dev.	.027 mg/L	
R.S.D.	12.5 %	

Precision of Duplicates-low range

	mid range	high range
s.d.	.015	1.004
mean	.670	1.106

W 0.01 mg/L

T 0.05 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- Extraction efficiency assumed to be 100% as material originally in solution. No ref. stds. available.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## BARIUM IN LANDFILL

Operating Range = 0.005 to 2.0 mg/L

### IN - RUN DUPLICATES

---

Range	<0.005	0.005 to 0.40	0.40 to 1.00	1.00 to 2.0	>2.0
no.	11	114	8	2	2
s.w.		0.0074	0.0145	1.0035	
mean		0.0852	0.6696	1.1060	

---

### QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	0.215	0.0268	12.50

---

### BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	84	0.005	0.0116

---

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Cadmium      TEST CODE: CDUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0

REVISION NO: Original

DATE:

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml

Container- Glass or plastic container

Preservative- do not preserve

Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-90%

Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.005 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.002 - 0.1 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

mean      1.760mg/L

std. dev.      0.129mg/L

R.S.D.      7.3 %

Precision of Duplicates-low range

mid range

high range

s.d.      0.0012

0.0003

mean      0.0068

0.0917

W .002 mg/L

T .010 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.

- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

CADMIUM

IN LANDFILL

Operating Range = 0.002 to 0.1 mg/L

## IN - RUN DUPLICATES

Range	<0.002	0.002 to 0.02	0.02 to 0.05	0.05 to 0.1	>0.1
no.	125	9	0	1	2
s.w.		0.0012	0.0000	0.0003	
mean		0.0068	0.0000	0.0917	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	1.760	0.1285	7.30

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	43	0.003	0.0036

DATE 87/05/15

5.10



ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Calcium                      TEST CODE: CAUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-100%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 50 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.02 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.2 - 100 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A

mean	304.1mg/L
std. dev.	34.2mg/L
R.S.D.	11.3 %

B

Precision of Duplicates-low range	mid range	high range
s.d.      0.4	0.6	2.2
mean      8.6	34.9	73.0

W 1.0 mg/L

T 5.0 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - assumed 100% as material in solution.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## CALCIUM IN LANDFILL

Operating Range = 0.200 to 100.0 mg/L

### IN - RUN DUPLICATES

Range	<0.200	0.200 to 20.00	20.00 to 50.00	50.00 to 100.0	>100.0
no.	3	14	16	45	59
s.w.		0.4014	0.5734	2.1665	
mean		8.5883	34.9423	72.9586	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	202	304.083	34.2314	11.26

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	43	0.279	0.3469

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Chromium      TEST CODE: CRUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-93%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.01 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

mean      1.73 mg/L  
std. dev.      0.15 mg/L  
R.S.D.      8.5 %

B

Precision of Duplicates-low range      mid range      high range

s.d.      0.012      0.025

mean      0.036      0.308

W 0.02 mg/L

T 0.10 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## CHROMIUM IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

### IN - RUN DUPLICATES

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	69	64	2	0	2
s.w.		0.0121	0.0253	0.0000	
mean		0.0359	0.3078	0.0000	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	1.731	0.1477	8.53

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	63	0.004	0.0044

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Cobalt                      TEST CODE: COUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-93%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.02 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A

mean	1.78 mg/L
std. dev.	0.10 mg/L
R.S.D.	5.6 %

B

Precision of Duplicates-low range                      mid range                      high range

s.d.                      0.007                      0.007

mean                      0.035                      0.405

W 0.02 mg/L

T 0.10 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.

- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

COBALT

IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

## IN - RUN DUPLICATES

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	111	24	1	0	1
s.w.		0.0067	0.0071	0.0000	
mean		0.0349	0.4050	0.0000	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	1.775	0.0994	5.60

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	47	0.002	0.0033

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Copper                      TEST CODE: CUUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AAO

REVISION NO: Original

DATE:

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-90%

Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.006 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.005 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A

B

mean                      1.84 mg/L

std. dev.                      0.11 mg/L

R.S.D.                      5.9 %

Precision of Duplicates-low range

mid range

high range

s.d.                      0.004

0.006

0.022

mean                      0.040

0.369

0.616

W .005 mg/L

T .025 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.

- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

COPPER

IN LANDFILL

Operating Range = 0.005 to 1.0 mg/L

## IN - RUN DUPLICATES

Range	<0.005	0.005 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	58	66	3	5	5
s.w.		0.0044	0.0059	0.0220	
mean		0.0401	0.3694	0.6161	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	1.838	0.1092	5.94

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	76	0.006	0.0094



ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Iron                      TEST CODE: FEUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-92%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 50 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.02 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.03 - 20.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A  
mean                      2.37 mg/L  
std. dev.                      0.26 mg/L  
R.S.D.                      11.0 %

B

Precision of Duplicates-low range	mid range	high range
s.d.                      0.11	0.05	0.32
mean                      0.60	5.94	15.04

W 0.1 mg/L

T 0.5 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference OCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.  
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

IRON

IN LANDFILL

Operating Range = 0.030 to 20.0 mg/L

## IN - RUN DUPLICATES

Range	<0.030	0.030 to 4.00	4.00 to 10.00	10.00 to 20.0	>20.0
no.	18	93	11	5	10
s.w.		0.1078	0.0510	0.3223	
mean		0.6019	5.9441	15.0433	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	197	2.370	0.2601	10.97

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	93	0.041	0.0516

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Lead                      TEST CODE: PBUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:  
Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-91%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.  
INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.02 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.03 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

mean                      .738 mg/L  
std. dev.                      .086 mg/L  
R.S.D.                      11.7 %

B

Precision of Duplicates-	low range	mid range	high range
s.d.	0.047		0.115
mean	0.078		0.594

W 0.05 mg/L

T 0.25 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

LEAD

IN LANDFILL

Operating Range = 0.030 to 1.0 mg/L

## IN - RUN DUPLICATES

Range	<0.030	0.030 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	121	11	0	2	3
s.w.		0.0472	0.0000	0.1149	
mean		0.0784	0.0000	0.5941	

## QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	201	0.738	0.0861	11.67

## BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	49	0.025	0.0329

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Magnesium      TEST CODE: MGUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original      DATE:  
NATURE OF LAST REVISION:

**SAMPLE HANDLING:**

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

**SAMPLE PREPARATION:** Partial Extn.-      Total Extn.-Yes % Extracted-100%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

**INTERFERENCES:** Several, compensated for by the computer program.

**REPORTING RESULTS:** mg/L to 3 places after decimal.

**INSTRUMENTATION:** Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 50 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.03 mg/L

**PERFORMANCE CHARACTERISTICS:**

Routine Operating Range- 0.2 - 100 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

	A	B
mean	288.0mg/L	
std. dev.	21.1mg/L	
R.S.D.	7.3 %	

Precision of Duplicates-	low range	mid range	high range
s.d.	0.30	0.66	4.00
mean	8.65	27.7	71.61

W 0.5 mg/L

T 2.5 mg/L

**CONTROL LIMITS:** Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

**REMARKS:**

- % extracted - assumed 100% as material in solution.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## MAGNESIUM IN LANDFILL

Operating Range = 0.200 to 100.0 mg/L

### IN - RUN DUPLICATES

---

Range	<0.200	0.200 to 20.00	20.00 to 50.00	50.00 to 100.0	>100.0
no.	6	64	25	21	21
s.w.		0.2951	0.6560	3.9991	
mean		8.6474	27.7172	71.6059	

---

### QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	288.016	21.1416	7.34

---

### BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	29	0.152	0.1474

---

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Manganese      TEST CODE: MNUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-94%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.003 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 2.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A  
mean      2.18 mg/L  
std. dev.      0.15 mg/L  
R.S.D.      6.6 %

B

Precision of Duplicates-low range

s.d.      0.005

mean      0.112

mid range

0.024

0.713

high range

0.033

1.366

W 0.01 mg/L

T 0.05 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference OCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.

- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## MANGANESE IN LANDFILL

Operating Range = 0.010 to 2.0 mg/L

### IN - RUN DUPLICATES

---

Range	<0.010	0.010 to 0.40	0.40 to 1.00	1.00 to 2.0	>2.0
no.	10	88	15	10	14
s.w.		0.0049	0.0240	0.0327	
mean		0.1115	0.7130	1.3657	

---

### QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	2.181	0.1447	6.64

---

### BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	48	0.004	0.0028

---



ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Molybdenum      TEST CODE: MOUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-100%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.005 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

mean	.677 mg/L
std. dev.	.108 mg/L
R.S.D.	15.9 %

B

Precision of Duplicates-low range	mid range	high range
s.d. 0.012		
mean 0.060		

W 0.01 mg/L

T 0.05 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF3 exceed  $\pm 15\%$  RSD.

REMARKS:

- Extraction efficiency assumed to be 100% as material originally in solution. No ref. stds. available.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## MOLYBDENUM IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

### IN - RUN DUPLICATES

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	123	13	0	1	0
s.w.		0.0120	0.0000	0.0243	
mean		0.0600	0.0000	0.5863	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	0.677	0.1077	15.91

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	87	0.012	0.0109

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Nickel                      TEST CODE: NIUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0  
REVISION NO: Original                      DATE:  
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-                      Total Extn.-Yes % Extracted-90%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.03 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

	A	B
mean	1.75 mg/L	
std. dev.	0.15 mg/L	
R.S.D.	8.6 %	

Precision of Duplicates-low range	mid range	high range
s.d. 0.006	0.015	0.012
mean 0.034	0.279	0.759

W 0.05 mg/L

T 0.25 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference OCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.

- Detection Limit - 3x std. dev. of low range within-run duplicates.  
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## NICKEL IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

### IN - RUN DUPLICATES

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	94	36	3	1	3
s.w.		0.0060	0.0148	0.0121	
mean		0.0338	0.2791	0.7587	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	204	1.752	0.1503	8.58

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	80	0.005	0.0056

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Vanadium      TEST CODE: VVUT      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0

REVISION NO: Original

DATE:

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-100%

Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.008 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.010 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

A

B

mean .587 mg/L

std. dev. .156 mg/L

R.S.D. 26.5 %

Precision of Duplicates-low range

mid range

high range

s.d. 0.059

0.015

mean 0.050

0.261

W 0.02 mg/L

T 0.10 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference QCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - assumed 100% as material in solution.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

## VANADIUM IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

### IN - RUN DUPLICATES

---

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	114	22	1	0	0
s.w.		0.0593	0.0151	0.0000	
mean		0.0503	0.2613	0.0000	

---

### QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	184	0.587	0.1557	26.54

---

### BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	48	0.009	0.0066

---

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: SELENIUM  
UNIT: Biomaterials

TEST CODE: SEUT

SAMPLE TYPE: Landfill Lch.  
SUPERVISOR: R. Sadana

METHOD CODE: 510BF3

REVISION NO:

DATE: January, 1983

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 100 ml

Container- Glass bottle with bakelite screw cap (16 oz.)

Preservative- 1ml conc. HNO<sub>3</sub> for sample filling 16 oz. bottle

Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-yes % Extracted-

Procedure- Digest samples by mixing 20 ml of sample and 4ml of acid (6:3:1 HNO<sub>3</sub>:HClO<sub>4</sub>:H<sub>2</sub>SO<sub>4</sub>) and heating on a hot plate until dense white fumes are formed. Allow sample to cool and add .5 ml of H<sub>2</sub>O and 2.5 ml HCl. Transfer the contents of the beaker to a test tube graduated at the 20 ml mark and bring to volume with washings. Shake and analyze by automated hydride-F.A.A.S. technique.

INTERFERENCES: Excessive concentrations of Cu, Fe, Ni

REPORTING RESULTS: mg/L-2 dec. if <10, 1 dec. if 10-100, 0 dec. if >100  
INSTRUMENTATION: Atomic Absorption Spectrophotometer (Varian 1200) with strip chart recorder, peristaltic pump, auto-sampler, open-ended and heated quartz 'T' cell (0.6x10 cm), and gas-liquid separator.

Calibration Range: 0 to 40 ng/ml

Resolution: 0.01 absorbance

Sensitivity: 20 ng/ml gives 0.200 Abs.

Instrument Detection Limit: 1 ng/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.001 to 0.04 mg/L

Accuracy-

Precision of Controls-

	A	B
mean	.075 mg/L	
std. dev.	.004 mg/L	
R.S.D.	5.3 %	

Precision of Duplicates-low range

mid range

high range

s.d. 0.0002

mean 0.002

W .001 mg/L

T .005 mg/L

CONTROL LIMITS:

REMARKS:

# SUMMARY REPORT OF QUALITY CONTROL DATA

## SELENIUM IN LANDFILL

Operating Range = 0.001 to 0.04 mg/L

### IN - RUN DUPLICATES

Range	<0.001	0.001 to 0.01	0.01 to 0.02	0.02 to 0.04	>0.04
no.	24	9	0	0	2
s.w.		0.0002	0.0000	0.0000	
mean		0.0020	0.0000	0.0000	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
475-3	34	0.075	0.0040	5.33

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
BLK	0	0.000	0.0000



ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Zinc                      TEST CODE: ZNUT                      SAMPLE TYPE: Landfill Lch.  
UNIT: Ind., Dom., Landfill Wastes                      SUPERVISOR: J. Pimenta

METHOD CODE: 533AA0

REVISION NO: Original

DATE:

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic container  
Preservative- do not preserve  
Other- Field filtered samples preferred

SAMPLE PREPARATION: Partial Extn.-      Total Extn.-Yes % Extracted-88%  
Procedure- Do not shake samples. Transfer 50 ml of sample to a test tube marked at 50 ml and 12.5 ml. Add 1 ml conc. HNO<sub>3</sub> and dry at 105°C in a forced air oven. Add 2.0 ml of aqua-regia to the test tube and digest for 2 hours. Make volume to 12.5 ml with distilled water and mix well. Determine the elemental concentration by ICP emission spectroscopy.

INTERFERENCES: Several, compensated for by the computer program.

REPORTING RESULTS: mg/L to 3 places after decimal.

INSTRUMENTATION: Inductively coupled plasma emission spectrometer, Jobin-Yvon JY48P, equipped with autosampler and DEC computer system for concentration print-out; Commodore Pet microcomputer interface to LIS.

Calibration Range: 0 - 10 mg/L

Resolution:

Sensitivity:

Instrument Detection Limit: 0.005 mg/L

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.01 - 1.0 mg/L

Accuracy- Not known; no reference standards available

Precision of Controls-

mean	1.75 mg/L
std. dev.	0.18 mg/L
R.S.D.	10.1 %

B

Precision of Duplicates-	low range	mid range	high range
s.d.	0.019	0.007	0.014
mean	0.045	0.272	0.920

W 0.02 mg/L

T 0.10 mg/L

CONTROL LIMITS: Analysis repeated if values for in-house reference OCLF1 exceed  $\pm 15\%$  RSD.

REMARKS:

- % extracted - based on spike recovery.
- Detection Limit - 3x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

# SUMMARY REPORT OF QUALITY CONTROL DATA

ZINC

IN LANDFILL

Operating Range = 0.010 to 1.0 mg/L

## IN - RUN DUPLICATES

---

Range	<0.010	0.010 to 0.20	0.20 to 0.50	0.50 to 1.0	>1.0
no.	52	65	7	1	12
s.w.		0.0186	0.0068	0.0141	
mean		0.0450	0.2715	0.9200	

---

## QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
QCLF-3	201	1.748	0.1769	10.12

---

## BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
rb	103	0.023	0.0245

---

ANALYTICAL PROCEDURE  
Inorganic Trace Contaminants Section

TEST NAME: Total cyanide TEST CODE: CCNAUR SAMPLE TYPE: Landfill Lch.  
UNIT: QC-Project SUPERVISOR: J. Hipfner

METHOD CODE: 001AC2

REVISION NO:

DATE

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- 500 ml  
Container- Glass or plastic (preferred)  
Preservative- NaOH  
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.- % Extracted-100

Procedure- The sample is first run to see if there is cyanide present. The is run directly by the automated high temperature distillation with 25% H<sub>3</sub>PO<sub>4</sub>-5% H<sub>3</sub>PO<sub>2</sub> followed by a colourimetric analysis with chloramine T -isonicotinic acid -barbituric acid method.

If the total cyanide is > .01 mg/L then 5 to 250 ml of sample is manually distilled with 30 ml of 15%(w/v) tartaric acid. The distillate is collected in 50 ml of 1N NaOH, and analyzed by the automated Technicon distillation system referred to above.

INTERFERENCES: SCN interference is removed by distillation.

Distillable organics may interfere; also S= at high levels.

REPORTING RESULTS: Mg/l CN: 3 decimal places up to 3 significant figs

INSTRUMENTATION: Technicon AAI continuous flow analyzer

including pump, colourimeter, appropriate autosampler and recorder.

High temperature distillation apparatus (Technicon). Manual dist. app

Calibration Range: 0 to 0.2 mg/l as CN

Resolution: 0.001 mg/l

Sensitivity:

Instrument Detection Limit: 0.001 mg/l

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.00100 to 0.200 mg/l

Accuracy- 100%

Precision of Controls-

	A	B
mean	.110 mg/L	0.059
std. dev.	.0027mg/L	0.0026
R.S.D.	2.45 %	4.41 %

Precision of Duplicates-low range mid range high range

s.d.

mean

W .001 mg/L

T .005 mg/L

CONTROL LIMITS:

REMARKS: Pure CN standards are recovered 100% during manual distillation. Complex cyanides can normally be expected to be recovered at 100%.

# SUMMARY REPORT OF QUALITY CONTROL DATA

## TOTAL CYANIDE IN LANDFILL

Operating Range = 0.001 to 0.4 mg/L

### IN - RUN DUPLICATES

Range	<0.001	0.001 to 0.08	0.08 to 0.2	0.2 to 0.4	>0.4
no.	0	2	0	0	0
s.w.		0	0	0	
mean		0.011	0	0	

### QA CONTROL SAMPLES

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
qc-a	146	0.149	0.0049	3.29
qc-b	146	0.018	0.0022	12.22

### BLANKS

BLANK I.D.	NO.	MEAN	STD. DEV.
BLK	146	0.001	0

DATE 88/06/01

## ANALYTICAL PROCEDURE

TEST NAME: Free cyanide  
UNIT: QC-Project

SAMPLE TYPE:Landfill Lch.

5.39

# SUMMARY REPORT OF QUALITY CONTROL DATA

## FREE CYANIDE IN LANDFILL

Operating Range = 0.001 to 0.4 mg/L

### IN - RUN DUPLICATES

---

Range	<0.001	0.001 to 0.08	0.08 to 0.2	0.2 to 0.4	>0.4
no.	0	2	0	0	0
s.w.		0	0	0	
mean		0.001	0	0	

---

### QA CONTROL SAMPLES

---

SAMPLE I.D.	NO.	MEAN	STD. DEV.	R.S.D.
qc-a	135	0.151	0.0062	4.11
qc-b	135	0.018	0.0022	12.22

---

### BLANKS

---

BLANK I.D.	NO.	MEAN	STD. DEV.
BLK	135	0.001	0

---



(6824)

MOE/ANN/ITC/ALRN

MOE/ANN/ITC/ALRN

Ontario Ministry of the En  
1986 annual quality

assurance alrn

c.1 a aa